



## **External Threaded Cryo Vials LD-ECV102L**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# External Threaded Cryo Vials LD-ECV102L

## Features:

- Equipped with Medical class silicon seal-ring
- Working temperature is -196 °C to 121 °C
- DNase, RNase and Pyrogenic free
- Pre-Sterilized
- Self-standing but Conical bottom
- Used for storing biological material, human or animal cells

## Specifications:

<b>Type</b>	Self-standing
<b>Capacity</b>	2.8 mL
<b>Packaging</b>	25 pcs/Bag
<b>Length of Tube</b>	61.50 mm
<b>Diameter of Tube</b>	12.50 mm
<b>Quantity per Case</b>	500