



## **Internal Threaded Cryo Vials LD-ICV105L**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Internal Threaded Cryo Vials LD-ICV105L

## Features:

- Working temperature is -196 °C to 121 °C
- Caps are made of Medical class Silicon
- DNase, RNase and Pyrogenic free
- Pre-Sterilized
- Self-standing but Conical bottom
- Used for storing biological material, human and animal cells

## Specifications:

<b>Type</b>	Self-standing
<b>Capacity</b>	5.0 mL
<b>Packaging</b>	25 pcs/Bag
<b>Length of Tube</b>	90.0 mm
<b>Diameter of Tube</b>	12.80 mm
<b>Quantity per Case</b>	500