



## **PCR Well Plates LD-PWP102L**

**Labodam Equipment Ltd.**

**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# PCR Well Plates LD-PWP102L

## Features:

- Maximum G-Force is 3500 RCF
- Material: Highest Grade Polypropylene
- Compatible with range of PCR and qPCR machines
- Thin-walled plates for optimum PCR and real-time (RT, qPCR) results
- Pre-sterilized
- RNase, DNase and pyrogenic free

## Specifications:

|                   |            |
|-------------------|------------|
| Type              | Non-Skirt  |
| Color             | Nature     |
| Capacity          | 0.2 mL     |
| Packaging         | 15 pcs/box |
| Quantity per Case | 150        |