



## **PCR Well Plates LD-PWP202L**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# PCR Well Plates LD-PWP202L

## Features:

- Maximum G-Force is 3500 RCF
- Material: Highest Grade Polypropylene
- Compatible with range of PCR and qPCR machines
- Thin-walled plates for optimum PCR and real-time (RT, qPCR) results
- Pre-sterilized
- RNase, DNase and pyrogenic free

## Specifications:

<b>Type</b>	Semi-Skirt
<b>Color</b>	Nature
<b>Capacity</b>	0.2 mL
<b>Packaging</b>	15 pcs/box
<b>Quantity per Case</b>	150