



## **Racked Bar-Coded Cryogenic Vials LD-BCV104L**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Racked Bar-Coded Cryogenic Vials LD-BCV104L

## Features:

- Material: Polypropylene
- Low Temperature Storage: -196 °C to 121 °C
- Sterilization Level: SAL 10<sup>-6</sup>
- Chemically resistant and scratch-proof 2D Laser Etched code on Bottom
- U-shaped for maximum working volume
- Barcoded rack for tracking
- Inner Diameter: 0.7 mm
- Cap Color: Transparent, White, Blue, Orange and Green
- DNase, RNase and Pyrogenic free

## Specifications:

<b>Height</b>	33.3 mm
<b>Barcode</b>	2D
<b>Cap Type</b>	Screw Cap
<b>Capacity</b>	1.4 mL
<b>Packaging</b>	96 pcs/Rack, 48 Racks/Case
<b>Working Volume</b>	1.15 mL
<b>Quantity per Case</b>	4608