



## **Fixed Volume Pipettes LD-FVP200L**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Fixed Volume Pipettes LD-FVP200L

## Features:

- Autoclavable at 121 °C
- Material: Polycarbonate
- Excellent durability and chemical resistance
- Ergonomic design provides excellent operating experience
- Easy-to-read volume display
- Trouble-free calibration and maintenance

## Specifications:

<b>Test Volume</b>	5 µL
<b>Accuracy error</b>	1.3 %
<b>Precision error</b>	1.2 %