



**Variable Volume Fully Autoclavable  
Pipettes LD-VVP106L**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Variable Volume Fully Autoclavable Pipettes LD-VVP106L

## Features:

- Autoclavable at 121 °C
- Material: Polycarbonate
- Available in 0.1 µL to 5000 µL volume
- Ergonomic design provides excellent operating experience
- Easy-to-read volume display
- Trouble-free calibration and maintenance

## Specifications:

<b>Increment</b>	1 µL
<b>Test Volume</b>	200 µL <hr> 100 µL <hr> 50 µL
<b>Volume Range</b>	50 - 200µL
<b>Accuracy error</b>	0.60 % <hr> 0.80 % <hr> 1.00 %
<b>Precision error</b>	0.15 % <hr> 0.30 % <hr> 0.40 %