



## **Variable Volume Pipette LD-VVP208L**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Variable Volume Pipette LD-VVP208L

## Features:

- Autoclavable at 121 °C
- Material: Polycarbonate
- Available in 0.1 µL to 10000 µL volume
- Ergonomic design provides excellent operating experience
- Easy-to-read volume display
- Trouble-free calibration and maintenance

## Specifications:

|                        |                                 |
|------------------------|---------------------------------|
| <b>Increment</b>       | 5 µL                            |
| <b>Test Volume</b>     | 1000 µL <hr> 500 µL <hr> 200 µL |
| <b>Volume Range</b>    | 200 - 1000µL                    |
| <b>Accuracy error</b>  | 0.60 % <hr> 0.70 % <hr> 0.90 %  |
| <b>Precision error</b> | 0.20 % <hr> 0.25 % <hr> 0.30 %  |