



## **Pipette Pump LD-PPP101L**

**Labodam Equipment Ltd.**

**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Pipette Pump LD-PPP101L

## Features:

- Material: Polypropylene
- Thumbwheel for precised aspiration or dispensing
- Resistant to acids and alkalis
- Available in three colors

## Specifications:

<b>Color</b>	Green
<b>Capacity</b>	10 mL
<b>Packaging</b>	Bulk Packing
<b>Quantity Per Case</b>	100