



Filter Pipette Tips LD-FPT102L

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Filter Pipette Tips LD-FPT102L

Features:

- Universal compatibility
- Material: Polypropylene
- Micro point design
- Pre sterilized
- RNase, DNase and Pyrogenic free
- Optional Color “ Yellow (FPT112L and FPT113L)
- Optional Color “ Blue (FPT114L, FPT115L, FPT116L and FPT117L)

Specifications:

Type	Short
Length	31.7 mm
Capacity	10 μ L
Packaging	96 Tips/Rack, 50 Racks/ Case
Quantity per Case	4800