



**Low Retention Filter Pipette Tips LD-  
LFT107L**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Low Retention Filter Pipette Tips LD-LFT107L

## Features:

- Universal compatibility
- Material: Polypropylene
- Micro point design
- Pre sterilized
- RNase, DNase and pyrogenic free
- Optional Color “ Yellow (LFT106L)
- Optional Color “ Blue (LFT107L and LFT108L)

## Specifications:

<b>Type</b>	Long
<b>Length</b>	86.4 mm
<b>Capacity</b>	1000 µL
<b>Packaging</b>	96 Tips/Rack, 50 Racks/ Case
<b>Quantity per Case</b>	4800