



**Ultra Low retention Pipette Tips LD-  
LRT113L**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Ultra Low retention Pipette Tips LD-LRT113L

## Features:

- Universal compatibility
- Material: Polypropylene
- Autoclavable at +121Â°C
- DNase, RNase and Human DNA free
- Micro Point Design
- Pre Sterilized
- Long reach design
- Optional Color â€“ Yellow (LRT106L, LRT107L, LRT108L, LRT109L, LRT110L and LRT111L)
- Optional Color â€“ Blue (LRT112L, LRT113L, LRT114L, LRT115L, LRT116L and LRT117L)

## Specifications:

<b>Type</b>	Long
<b>Capacity</b>	1000 ?L
<b>Packaging</b>	1000 Tips/Bag, 10 Bags/Case
<b>Length of tube</b>	86.40 mm
<b>Quantity per Case</b>	10000