



## **Pipette Tip for Eppendorf LD-SPT105L**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Pipette Tip for Eppendorf LD-SPT105L

## Features:

- Material: Polypropylene
- Autoclavable
- Smooth inner wall to avoid leakage and sample residual
- Pre Sterilized
- DNase, RNase and pyrogen free

## Specifications:

<b>Type</b>	Filter
<b>Color</b>	Blue
<b>Capacity</b>	1000 $\mu$ L
<b>Packaging</b>	500 pcs/Pack, 20 Packs/Case
<b>Compatibility</b>	Eppendorf
<b>Quantity Per Case</b>	10000