



## **Universal Pipette Tips LD-UPT116L**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Universal Pipette Tips LD-UPT116L

## Features:

- Universal compatibility
- Material: Polypropylene
- Autoclavable at +121Â°C
- DNase, RNase and Human DNA free
- Micro Point Design
- Pre Sterilized
- Long reach design
- Optional Color â€“ Yellow (UPT106L, UPT107L, UPT108L, UPT109L, UPT110L and UPT111L)
- Optional Color â€“ Blue (UPT112L, UPT113L, UPT114L, UPT115L, UPT116L, UPT117L and UPT118L)

## Specifications:

<b>Type</b>	Long
<b>Capacity</b>	1250/1500 ?L
<b>Packaging</b>	500 Tips/Bag , 20 Bags/Case
<b>Length of tube</b>	103.20 mm
<b>Quantity per Case</b>	10000