



## **Polystyrene Serological Pipette LD- PSP100L**

**Labodam Equipment Ltd.**  
**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# Polystyrene Serological Pipette LD-PSP100L

## Features:

- Good optical Transmittance
- Material: Polystyrene
- Gamma ray sterilized
- Equipped with Porex filter
- DNase, RNase and endotoxin free

## Specifications:

<b>Capacity</b>	1 mL
<b>Packaging</b>	1000 Pcs/Case
<b>Color of Mark</b>	Yellow
<b>Min. graduation</b>	1/100 mL
<b>Quantity per Case</b>	1000