



Polystyrene Serological Pipette LD- PSP102L

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Polystyrene Serological Pipette LD-PSP102L

Features:

- Good optical Transmittance
- Material: Polystyrene
- Gamma ray sterilized
- Equipped with Porex filter
- DNase, RNase and endotoxin free

Specifications:

Capacity	5 mL
Packaging	500 Pcs/Case
Color of Mark	Blue
Min. graduation	1/10 mL
Quantity per Case	500