



Polystyrene Serological Pipette LD- PSP105L

Labodam Equipment Ltd.
info@labodam.com | www.labodam.com

Polystyrene Serological Pipette LD-PSP105L

Features:

- Good optical Transmittance
- Material: Polystyrene
- Gamma ray sterilized
- Equipped with Porex filter
- DNase, RNase and endotoxin free

Specifications:

Capacity	50 mL
Packaging	100 Pcs/Case
Color of Mark	Black
Min. graduation	5/10 mL
Quantity per Case	100