



## **PolyEthylene Transfer Pipette LD-PTP100L**

**Labodam Equipment Ltd.**

**[info@labodam.com](mailto:info@labodam.com) | [www.labodam.com](http://www.labodam.com)**

# PolyEthylene Transfer Pipette LD-PTP100L

## Features:

- Material: Low-density Polyethylene
- Heat seal on the tip can prevent liquid leakage
- Gamma ray sterilized
- Non-toxic
- High Transparency and Elasticity

## Specifications:

<b>Length</b>	143 mm
<b>Capacity</b>	1 ml
<b>Drop Vol.</b>	42 ?L
<b>Packaging</b>	500 pcs/pk
<b>Quantity Per Case</b>	5000