



114 L Flammable Storage Cabinet
LFSC-A10

Description :

114 L Flammable Storage Cabinet LFSC-A10 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Specifications :

Capacity	114 L
Door type	Dual door manual
Adjustable shelf	1
Shelves loading capacity	50 kgs
Dimension (H x W x D)	1120 x 1090 x 460 mm
Weight	107 kgs

Features :

- Capacity - 114 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- One adjustable shelf
- Adjustable spill retaining shelves

Applications :

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)



Labodam Equipment Ltd

18a Melton Road Leicester LE4 5EA
United Kingdom

www.labodam.com // info@labodam.com