



**227 L COMBUSTIBLE CABINET LD-LCBC-A12**

**[www.labodam.com](http://www.labodam.com) | [info@labodam.com](mailto:info@labodam.com)**

227 L Combustible Cabinet LD-LCBC-A12 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three-point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead-free powder coating provides durability and chemical resistance preventing corrosion and humidity.

## Features

- Capacity - 227 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three-point latch and manual lock
- Two adjustable shelves
- Adjustable spill retaining shelves

## Applications

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

## Specifications

Model No.	LD-LCBC-A12
Capacity	227 L
Door type	Dual door manual
Adjustable shelf	2
Shelves loading capacity	100 kgs
Dimension (H x W x D)	1650 x 860 x 860 mm
Weight	100 kgs