



45 L Flammable Storage Cabinet
LFSC-B12

Description :

45 L Flammable Storage Cabinet LFSC-B12 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

Specifications :

Capacity	45 L
Door type	Single door manual
Adjustable shelf	1
Shelves loading capacity	50 kgs
Dimension (H x W x D)	890 x 590 x 460 mm
Weight	59 kgs

Features :

Capacity - 45 L

Dual walled manual

38 mm insulating air space

180° fully opening door

Three point latch and manual lock

One adjustable shelf

Adjustable spill retaining shelves

Applications :

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)



Labodam Equipment Ltd

18a Melton Road Leicester LE4 5EA
United Kingdom

www.labodam.com // info@labodam.com