

BLOOD BAG TUBE SEALER

LD-LBTS-A10



BLOOD BAG TUBE SEALER LBTS-A10

Blood Bag Tube Sealer LD-LBTS-A10 is compact unit used to seal the blood bags, infusion bags and urine bags with high frequency sealing system. It is suitable for blood bags with diameter 3 to 6 mm and 0.5 to 2 s sealing time. Designed with single chip and automatic sealing voltage adjusted. Equipped with electronic tube of about 1000 hours lifespan.

Features:

- Designed with single chip and automatic sealing voltage adjusted
- Used in blood collecting vehicle with suction feet to avoid moving
- Sealing time ranges from 0.5 to 2 s
- 1000 hours lifespan of electronic tube
- Instantaneous sealing with uniform heating, firm sealing, pollution-freely
- Designed with circulating air duct, effectively reducing the inside temperature and constantly work
- Sealer maintains quality of key components avoids hemolysis of blood and ensure high quality

Applications:

Blood Bag Tube Sealer is used in blood bags, infusion bags and urine bags by high frequency sealing system. Also used in blood-bank, hospital, biological industry, etc.

Specifications:

Model No.	LD-LBTS-A10
Sealing Time	0.5 to 2 s
Sealing Type	Automatic
Diameter of tube	Φ3 to 6 mm
Lifespan of Electronic Tube	About 1000 hours
External dimension	340 × 180 × 170 mm
Packaging dimension	400 × 280 × 270 mm
Net weight	9 kg
Gross weight	12 kg