



Conical Bottom Centrifuge Tubes
CCT100L

Description :

Specifications :

Capacity	15 mL
Length of tube	12 cm
Diameter outside tube	17 mm
Diameter of Cap	22 mm
Packaging	25 pcs/Bag
Quantity per Case	500

Features :

Maximum RCF 12000 x g

Material: Virgin Polypropylene (PP)

Working Temperature is -80 °C to 121 °C

Transparent wall permits easy observation of tube contents

Graduation in white ink

Frosted marking area

Flat-top screw cap for free handling



Labodam Equipment Ltd

18a Melton Road Leicester LE4 5EA

United Kingdom

www.labodam.com // info@labodam.com