



Conical Bottom Centrifuge Tubes
CCT101L

Description :

Specifications :

Capacity	15 mL
Length of tube	12 cm
Diameter outside tube	17 mm
Diameter of Cap	22 mm
Packaging	50 pcs/Foam Rack
Quantity per Case	500

Features :

- Maximum RCF 12000 x g
- Material: Virgin Polypropylene (PP)
- Working Temperature is -80 Â°C to 121 Â°C
- Transparent wall permits easy observation of tube contents
- Graduation in white ink
- Frosted marking area
- Flat-top screw cap for free handling



Labodam Equipment Ltd

18a Melton Road Leicester LE4 5EA
United Kingdom

www.labodam.com // info@labodam.com