



*Conical Bottom Centrifuge Tubes*  
*CCT200L*

## Description :

## Specifications :

Capacity	50 mL
Length of tube	11.7 cm
Diameter outside tube	30 mm
Diameter of Cap	34 mm
Packaging	25 pcs/Bag
Quantity per Case	500

## Features :

Maximum RCF 12000 x g

Material: Virgin Polypropylene (PP)

Working Temperature is -80 Â°C to 121 Â°C

Transparent wall permits easy observation of tube contents

Graduation in white ink

Frosted marking area

Flat-top screw cap for free handling

The logo for Labodam features the word "labodam" in a bold, black, sans-serif font. The letter "a" is stylized with a green circular graphic element inside it, consisting of a solid green circle at the bottom and a dotted green circle at the top.

**Labodam Equipment Ltd**

18a Melton Road Leicester LE4 5EA

United Kingdom

[www.labodam.com](http://www.labodam.com) // [info@labodam.com](mailto:info@labodam.com)