



*Filter Pipette Tips FPT104L*

## Description :

## Specifications :

Capacity	20 µL
Length	50 mm
Packaging	96 Tips/Rack, 50 Racks/ Case
Quantity per Case	4800

## Features :

Universal compatibility

Material: Polypropylene

Micro point design

Pre sterilized

RNase, DNase and Pyrogenic free

Optional Color “ Yellow (FPT112L and FPT113L)

Optional Color “ Blue (FPT114L, FPT115L, FPT116L and FPT117L)



**Labodam Equipment Ltd**

18a Melton Road Leicester LE4 5EA  
United Kingdom

[www.labodam.com](http://www.labodam.com) // [info@labodam.com](mailto:info@labodam.com)