



Filter Pipette Tips FPT112L

www.labodam.com // info@labodam.com

Description:

Specifications:

Length	58.7 mm
Capacity	200 µL
Packaging	96 Tips/Rack, 50 Racks/ Case
Quantity per Case	4800

Features:

- Universal compatibility
- Material: Polypropylene
- Micro point design
- Pre sterilized
- RNase, DNase and Pyrogenic free
- Optional Color “ Yellow (FPT112L and FPT113L)
- Optional Color “ Blue (FPT114L, FPT115L, FPT116L and FPT117L)



Labodam Equipment Ltd

18a Melton Road Leicester LE4 5EA

United Kingdom

www.labodam.com // info@labodam.com