



*Fixed Volume Fully Autoclavable
Pipettes FVP101L*

Description :

Specifications :

Test Volume	10 μ L
Accuracy error	0.8 %
Precision error	0.8 %

Features :

Autoclavable at 121 $^{\circ}$ C

Material: Polycarbonate

Excellent durability and chemical resistance

Ergonomic design provides excellent operating experience

Easy-to-read volume display

Trouble-free calibration and maintenance



Labodam Equipment Ltd

18a Melton Road Leicester LE4 5EA

United Kingdom

www.labodam.com // info@labodam.com