



*Fixed Volume Fully Autoclavable
Pipettes FVP105L*

Description :

Specifications :

Test Volume	100 µL
Accuracy error	0.5 %
Precision error	0.3 %

Features :

Autoclavable at 121 Â°C

Material: Polycarbonate

Excellent durability and chemical resistance

Ergonomic design provides excellent operating experience

Easy-to-read volume display

Trouble-free calibration and maintenance

The logo for Labodam features a stylized green icon of a laboratory flask or beaker with three dots above it, followed by the word "labodam" in a bold, black, sans-serif font.

Labodam Equipment Ltd

18a Melton Road Leicester LE4 5EA

United Kingdom

www.labodam.com // info@labodam.com