



Lab Refrigerator-Freezer Combination LRFC-A11

www.labodam.com // info@labodam.com

Description :

Lab Refrigerator-Freezer Combination LRFC-A11 is an upright type model designed for laboratory use offering superior product preservation with exceptional product quality. The four units caster facilitate easy handling. It also has various safety features like high or low temperature alarm, system failure audible/visual alarm. The temperature ranges from 2 °C to 8 °C for refrigerator and (-10) °C to (-25) °C for freezer. Also has wide range of capacity as 225 L for refrigerator and 225 L for freezer.

Specifications :

Effective capacity	450 Liters
Refrigerator room	225 Liters
Freezer room	225 Liters
Temperature range: Up room:	2 °C to 8 °C
Temperature range: Down room:	(-10) °C to (-25) °C
External size (W x D x H)	670 x 780 x 1850 mm
Inside size (W x D x H)	640 x 485 x 655 mm : UC 640x530x655 mm : LC
Input power	320 W
Power supply	220 V / 50 Hz, 220 V 60 Hz or 110 V 50 / 60 Hz
Net Weight	125 kg
Gross Weight	140 kg

Features :

Microprocessor temperature controller
LCD display
Upright type
Interior made up of ABS material
Painted steel board outer body
Safety lock on the door with door handle for additional security
Easy handling facilitated by 4 units caster
7-12units drawers made of plastic inside
R507, CFC free refrigerant
Rapid cool-down and very fast recovery due to fan assistance
High or low temperature alarm, system failure audible/visual alarm
Optional: chart recorder

Applications :

Laboratory refrigerators and freezers are generally used for storing non-volatile reagents and non-volatile biological specimens and used for storing liquid drugs, vaccines, samples in laboratories and research centres.



Labodam Equipment Ltd

18a Melton Road Leicester LE4 5EA

United Kingdom

www.labodam.com // info@labodam.com