



*LOW RETENTION FILTER PIPETTE TIPS
LD-LFT103L*

www.labodam.com | info@labodam.com

Specifications:

Model No.	LD-LFT103L
Length	50 mm
Capacity	50 µL
Packaging	96 Tips/Rack, 50 Racks/ Case
Quantity per Case	4800

Features:

- Universal compatibility
- Material: Polypropylene
- Micro point design
- Pre-sterilized
- RNase, DNase and pyrogenic free
- Optional Color Yellow (LFT106L)
- Optional Color Blue (LFT107L and LFT108L)


Labodam Equipment Ltd.

18A Melton Road, Leicester,
LE4 5EA, United Kingdom

www.labodam.com | info@labodam.com