

MAGNETIC STIRRERS

LMAS-A1 Series



Magnetic Stirrer LMAS-A10

Magnetic Stirrer LMAS-A10 is a microprocessor controlled stirrer equipped with a stirring capacity of 3 L and a speed range of 150 to 1500 rpm. Equipped with LED display for convenient control and monitoring of the stirring speed. Features individual control for each plate with an easy to clean anti-corrosive stainless steel plate surface.

Features

- Microprocessor controlled unit
 - LED display (control and monitoring of stirring speed)
 - Stainless steel stirring plate surface
 - Individual control for each plate
-

Magnetic Stirrer LMAS-A11

Magnetic Stirrer LMAS-A11 is a microprocessor controlled stirrer equipped with a stirring capacity of 2 L and a speed range of 150 to 1500 rpm. Equipped with LED display for convenient control and monitoring of the stirring speed. Features individual control for each plate with an easy to clean anti-corrosive stainless steel plate surface.

Features

- Microprocessor controlled unit
 - LED display (control and monitoring of stirring speed)
 - Stainless steel stirring plate surface
 - Individual control for each plate
-

Magnetic Stirrer LMAS-A12

Magnetic Stirrer LMAS-A12 is a microprocessor controlled stirrer equipped with a stirring capacity of 2 L and a speed range of 150 to 1500 rpm. Equipped with LED display for convenient control and monitoring of the stirring speed. Features individual control for each plate with an easy to clean anti-corrosive stainless steel plate surface.

Features

- Microprocessor controlled unit
- LED display (control and monitoring of stirring speed)
- Stainless steel stirring plate surface
- Individual control for each plate

Applications

Used to mix component, either solids or liquids and get a homogeneous liquid mixture etc.

Specifications

Model No	LMAS-A10	LMAS-A11	LMAS-A12
Mixing capacity	3 L	2 L	2 L
Number of plates	1	4	6
Plate dimension	158 x 143 mm	313 x 303 mm	473 x 303 mm
Speed range	150 rpm ~ 1500 rpm	150 rpm ~ 1500 rpm	150 rpm ~ 1500 rpm
Display mode	LED	LED	LED
Overall dimensions	189 x 215 x 48 mm	344 x 377 x 48 mm	504 x 377 x 48 mm