



Manual Dual Door  
**Corrosive Cabinet**  
LD-LCOC-A1 SERIES



### 114 L Corrosive Cabinet LD-LCOC-A10

114 L Corrosive Cabinet LD-LCOC-A10 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features one adjustable shelf to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

#### Features

- Capacity - 114 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- One adjustable shelf
- Adjustable spill retaining shelves

### 114 L Corrosive Cabinet LD-LCOC-A11

114 L Corrosive Cabinet LD-LCOC-A11 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

#### Features

- Capacity - 170 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- Two adjustable shelves
- Adjustable spill retaining shelves

#### Applications

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

## Specifications

Model No	LD-LCOC-A10	LD-LCOC-A11
Capacity	114 L	170 L
Door type	Dual door manual	Dual door manual
Adjustable shelves	1	2
Shelves loading capacity	50 kgs	100 kgs
Dimension ( H x W x D )	1120 x 1090 x 460 mm	1650 x 1090 x 460 mm
Weight	109 kgs	151 kgs

### 114 L Corrosive Cabinet LD-LCOC-A12

114 L Corrosive Cabinet LD-LCOC-A12 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

### Features

- Capacity - 227 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- Two adjustable shelves
- Adjustable spill retaining shelves

### 114 L Corrosive Cabinet LD-LCOC-A13

114 L Corrosive Cabinet LD-LCOC-A13 is a dual walled manual fire resistance cabinet with 38 mm of insulating air space. Doors are reinforced with a three point latch and manual lock for security. Features two adjustable shelves to accommodate storage for different containers with 2 inches vents integrated with flame arresters. Equipped with spill catcher shelves for preventing incidental drips. Inner and outer lead free powder coating provides durability and chemical resistance preventing corrosion and humidity.

#### Features

- Capacity - 340 L
- Dual walled manual
- 38 mm insulating air space
- 180° fully opening door
- Three point latch and manual lock
- Two adjustable shelves
- Adjustable spill retaining shelves

#### Applications

Used for storage of flammable and combustible liquids away from incompatible materials (e.g., oxidizers)

#### Specifications

Model No	LD-LCOC-A12	LD-LCOC-A13
Capacity	227 L	340 L
Door type	Dual door manual	Dual door manual
Adjustable shelves	2	2
Shelves loading capacity	100 kgs	100 kgs
Dimension ( H x W x D )	1650 x 860 x 460 mm	1650 x 1090 x 860 mm
Weight	151 kgs	186 kgs