



# MEDICAL REFRIGERATOR LD-LMR-A10

# MEDICAL REFRIGERATOR LD-LMR-A10

The Medical Refrigerator LD-LMR-A10 is an upright type freezer comes with 260 L capacity. Offers 2°C to 8°C adjustable temperature range with increment at 0.1°C and digital temperature displays with 0.1°C resolution Incorporated with microprocessor controlled system with interior LED display for ease of monitoring. Equipped with eco- friendly refrigerant HCFC-Free and CFC-freezer. Designed with forced air cooling system and auto defrost. Built in with exterior epoxy coated, anti-bacteria and easy cleaning. With 4 shelves design and 4 casters which aids in easy mobility and handling The spirit pharmacy refrigerator is suitable for installation in drug stores, pharmacy companies, hospitals, epidemic prevention centres to store medicines, testing kits, reagents, vaccines etc.

## FEATURES

- ◆ Designed with High-tech integrated sensors to display and control temperature
- ◆ Single glass door design with safety lock for safe operation
- ◆ Removable outdoor double seal design for better temperature maintain and energy saving
- ◆ Forced air cooling system
- ◆ Lockable caster for easy mobility
- ◆ Condensation on glass door is prevented by electrical heater
- ◆ Shelves adjustable design
- ◆ Password protection to prevent unauthorized temperature setting
- ◆ Built- in alarm for High and low temperature, power failure, sensor error alarm, low battery alarm, door ajar alarm
- ◆ Incorporated with buzzer and flash alarm styles and start delay protection
- ◆ It is complied with Certificate CE, ISO9001

## APPLICATIONS

The Medical Refrigerator is designed to store reagent, vaccines, blood plasma, and many other biological materials which is widely used in research, electronic, clinics, institution and medical fields. It is suitable for use in culture media, reagents and samples for testing, as well as for storage of enzymes for genetic research.

## SPECIFICATION

<b>Model No.</b>	<b>LD-LMR-A10</b>
Temperature Range	2 to 8°C
Uniformity	± 3°C
Cabinet Type	Upright
Gross Capacity	260 L
Climate Class	10 to 32°C
Cooling Type	Forced Air Cooling
Defrost Mode	Automatic
Refrigerant	CFC-Free
Door	Single Glass Door
Electrical Heater	Optional
Noise	42 dB
Temperature Control	Microprocessor
Display	LED
Self	4
Power Supply	220 V to 240 V/50 Hz to 60 Hz
Power	254 W
Internal Dimension (W×D×H)	500×420×1090 mm
External Dimension (W×D×H)	580×560×1785 mm
Packing Dimension (W×D×H)	710×665×1900
Net Weight	63 kg
Gross Weight	65 kg

## OPTIONAL ACCESSORIES

<b>Sr.No</b>	<b>Name of Accessories</b>	<b>Quantity</b>
1.	Test Hole	1
2.	Printer	1
3.	Chart Recorder	1
4.	USBI Interface	1
5.	Data Logger	1