



## MEDICAL SUCTION MACHINE

# LD-LMSM-A10

## Medical Suction Machine LD-LMSM-A10

Medical suction machine is designed for wound care with constant and random mode. It is applicable in skin care and wound management. The suction filter prevents any foreign contaminants from pump. Analogue pressure gauge for indication of negative pressure value. Its portable feature and lightweight design can be operated in operating room and also helpful in other necessary locations.

### Features

- Portable design convenient in operation in emergency room and other remote location
- Two working modes: Continuity and Intermittence can switch between modes
- Liquid collection vessel made of polycarbonate for collection of fluids up to 1200 mL
- Suction filter for prevention of any contaminants
- Manual knob for precise vacuum control
- Analogue gauge for indication of negative pressure in mmHg and MPa

### Application

Medical Suction Machine is used in emergency room for wound management and skin care to remove excess fluids, secretions, or foreign materials from a patient's airway or body cavity. They play a crucial role in various medical settings, from hospitals and clinics to emergency medical services and home care

### Specifications

| Model No.                       | LD-LMSM-A10                              |
|---------------------------------|------------------------------------------|
| Mode                            | Continuity / Intermittence               |
| Suction Pump                    | Piston                                   |
| Maximum Negative Pressure Value | $\leq 0.06$ MPa                          |
| Negative Pressure               | 0.02 MPa - 0.06 MPa                      |
| Pumping Rate                    | $\leq 20$ L/min                          |
| Time Adjustment                 | Working time adjustment range: 1-10 min  |
|                                 | Interval time adjustment range: 1-10 min |
|                                 | Default time: 5 min (ON)/2 min (OFF)     |
| Container                       | 1200 mL                                  |

**Medical Suction Machine LD-LMSM-A10**

|              |                  |
|--------------|------------------|
| Noise        | ≤ 55dB (A)       |
| Power Supply | AC 220V 50-60 Hz |
| Input Power  | 100 VA           |