



*Microcentrifuge Tubes MCT300L*

## Description :

## Specifications :

Capacity	2.0 mL
Type	Common
Length of tube	42 mm
Diameter of the tube	11 mm
Packaging	500 pcs/Bag
Quantity per Case	5000

## Features :

Maximum RCF 28000 x g

Working Temperature is -90 °C to 121 °C

Conical bottom design

Autoclavable temperature is +121 °C

Material: Polypropylene

Frosted marking area on tube and cap

Attached cap



**Labodam Equipment Ltd**

18a Melton Road Leicester LE4 5EA  
United Kingdom

[www.labodam.com](http://www.labodam.com) // [info@labodam.com](mailto:info@labodam.com)