



PCR Tube PTB100L

Description :

Specifications :

Specification	0.2 mL
Cap Type	Flat Cap
Packaging	1000 pcs/Bag, 70 Bags/Case
Quantity per case	70000

Features :

Material: Polypropylene

Ultra-thin walls for efficient thermal transfer and maximum yield

Domed or flat insertion caps ensure a tight fit that prevents sample evaporation

Optically clear cap for real time qPCR



Labodam Equipment Ltd

18a Melton Road Leicester LE4 5EA

United Kingdom

www.labodam.com // info@labodam.com