



PCR Well Plates PWP103L

Description :

Specifications :

Capacity	0.2 mL
Type	Non-Skirt
Color	White
Packaging	15 pcs/box
Quantity per Case	150

Features :

Maximum G-Force is 3500 RCF

Material: Highest Grade Polypropylene

Compatible with range of PCR and qPCR machines

Thin-walled plates for optimum PCR and real-time (RT, qPCR) results

Pre-sterilized

RNase, DNase and pyrogenic free



Labodam Equipment Ltd

18a Melton Road Leicester LE4 5EA
United Kingdom

www.labodam.com // info@labodam.com