



Plastic Pipette Stand PPS100L

Description :

Specifications :

| | |
|-------------------|--------------------|
| Type | Z Type |
| Material | PMMA |
| Holds | 5 |
| Color | White / Blue / Red |
| Packaging | Bulk Packing |
| Quantity Per Case | 10 |

Features :

Material: Polymethyl-methacrylate and Acrylonitrile butadiene styrene

Available in 3 different models

Ensures easy removal and replacement of pipettes with weighted base

Convenient and safe storage

The logo for Labodam features the word "labodam" in a bold, black, sans-serif font. The letter "a" is replaced by a stylized green icon of a laboratory flask with three small green dots above it, representing a pipette tip or a small amount of liquid.

Labodam Equipment Ltd

18a Melton Road Leicester LE4 5EA

United Kingdom

www.labodam.com // info@labodam.com