



*Polystyrene Serological Pipette*  
*PSP103L*

## Description :

## Specifications :

Capacity	10 mL
Min. graduation	1/10 mL
Color of Mark	Orange
Packaging	500 Pcs/Case
Quantity per Case	500

## Features :

Good optical Transmittance

Material: Polystyrene

Gamma ray sterilized.

Equipped with Porex filter

DNase, RNase and endotoxin free



**Labodam Equipment Ltd**

**18a Melton Road Leicester LE4 5EA**

**United Kingdom**

**[www.labodam.com](http://www.labodam.com) // [info@labodam.com](mailto:info@labodam.com)**