



*Screw Cap Microtubes SCM106L*

## Description :

## Specifications :

|                   |                           |
|-------------------|---------------------------|
| Capacity          | 2.0 mL                    |
| Cap Type          | O-ring                    |
| Cap Assembled     | No                        |
| Packaging         | 500 Pcs/Bag, 20 Bags/Case |
| Quantity per Case | 10000                     |

## Features :

Material: Polypropylene

Equipped with O-ring for good sealability

Conical Bottom

Pre Sterilized

Free of DNase and RNase

Available in Yellow, Blue, Pink, Purple and Amber color



**Labodam Equipment Ltd**

18a Melton Road Leicester LE4 5EA  
United Kingdom

[www.labodam.com](http://www.labodam.com) // [info@labodam.com](mailto:info@labodam.com)