



*Single Channel Electronic pipette*  
**SEP100L**

## Description :

## Specifications :

Volume Range	0.5 - 10 $\mu\text{L}$
Increment	0.01 $\mu\text{L}$
Test Volume	10 $\mu\text{L}$
	1 $\mu\text{L}$
Accuracy Error	$\pm 1.00\%$
	$\pm 3.50\%$
Precision Error	0.50 %
	3.00 %

## Features :

Dual charging modes for longer operation time

Digital control and multifunctioning Pipette

Autoclavable (Lower part)

Adjustable speeds for aspiration and dispensing

Self-calibration, can be connected to a PC via USB for calibration

Accuracy Error and Precision Error is given according to ISO8655-2 standard



**Labodam Equipment Ltd**

**18a Melton Road Leicester LE4 5EA**

**United Kingdom**

**[www.labodam.com](http://www.labodam.com) // [info@labodam.com](mailto:info@labodam.com)**