



Single Punch Tablet Press Machine LSTP-A1 Series

Single Punch Tablet Press Machine LSTP-A1 Series

Single Punch Tablet Press Machine LSTP-A1 Series features small desktop type press for motive and control sheeting. Mounted with one set of punch and die, the depth of filling material tablet press and thickness of the tablet is adjustable. Tablets of different granular materials are pressed under pressure for use in research and quality preparations across pharmaceutical, food, chemical industry, cosmetic industry etc.

Features:

- ❏ Desktop type press for motive and control sheeting
- ❏ Adjustable depth of filling material for tablet pressing
- ❏ Adjustable thickness of the tablet
- ❏ Minimum consumption of lab material 200 g
- ❏ High precision of filling
- ❏ Low noise operations.
- ❏ Hand and motor operations (LSTP-A11, LSTP-A12) , Only motor operations (LSTP-A10)
- ❏ Easy to use and maintain

Application:

Used in pharmaceutical preparations, food, chemical and other industries for research and quality analysis.

Specifications

Model no.	LSTP-A10	LSTP-A11	LSTP-A12
Maximum pressure	60 KN	50 KN	15 KN
Maximum diameter Tablet	25 mm	20 mm	12 mm
Maximum depth of fill	16 mm	12 mm	11 mm
Maximum thickness of tablet	10 mm	6 mm	
Production capacity	4000 pcs/ hour	5000 pcs/ hour	
Dies	1 set of round shape dies		
Motor power	1100 W	750 W	550 W
Motor speed	1400 rpm		
Net weight	160 kg	125 kg	85 kg
HS Code	8462991000		