



TABLE TOP AUTOCLAVE

LT TA-A20

Table top autoclave LTTA-A20

Table top autoclave LTTA-A20 is a robust and compact, space efficient unit with 20 L of capacity and a sterilizing temperature of 134 °C. It is equipped with a key feature of 4 - 6 mins rapid sterilization along with safety features like automatic shutoff preventing overheating of chamber, automatic cold air discharge and auto protection against raised temperature and pressure. Integrated with steam-water inner circulation system with steam discharge making the environment for sterilization clean and dry

Features

- Digital LCD display to monitor working status and parameters
- Working temperature : 134 °C
- 4 - 6 mins rapid sterilization
- Automated drying function
- Equipped with interior circulate system
- Inner steam circulatory system
- Automatic cold air discharge from system chamber ensures dry and clean environment
- Door lock feature : Door cannot be opened till the pressure inside the chamber doesn't fall to 0.027 MPa
- Self-inflating leak proof chamber
- Auto exhaustion of cool air
- Safe door lock system ensuring safety / protection during sterilization
- Automatic power cut off with alarm with low water levels indication
- Stainless steel sterilizing baskets

Application

Used to sterilize culture media, instruments, dressings, intravenous equipments, rubber closure of glass equipments used in aseptic techniques and numerous objects that can withstand high temperature and pressure

Specification

Model	LTTA-A20
Capacity	24 L
Number of sterilizing trays	3
Sterilization Temperature	134 °C
Sterilization Pressure	0.22 Mpa
Timer Range	0-60 Minutes
Heat Average	$\leq \pm 1$ °C
Chamber Dimension ($\Phi \times D$)	270×410 mm
Net Weight	18 Kgs
Gross Weight	20 Kgs
Power Supply	AC 220 V 50 Hz/1.5 KW