



# TABLE TOP AUTOCLAVE

## LTTA-A1 SERIES

## Table top Autoclave LTTA-A1 Series

Table top autoclave LTTA-A1 series is robust and compact, space efficient unit with a sterilizing temperature of 134°C. It is equipped with 3 sterilizing trays along with safety features like automatic shutoff preventing overheating of chamber, automatic cold air discharge and auto protection against raised temperature and pressure making a reliable and user friendly device.

### Features

- ▶ Working temperature : 134 °C
- ▶ Door lock feature : Door cannot be opened till the pressure inside the chamber doesn't fall to 0.027 MPa
- ▶ Safe door lock system ensuring safety / protection during sterilization
- ▶ Automatic power cut off with alarm with low water levels indication
- ▶ Stainless steel sterilizing trays
- ▶ Automatic cold air discharge from system chamber ensures dry and clean environment
- ▶ Auto exhaustion of cool air

### Application

Used to sterilize culture media, instruments, dressings, intravenous equipments, rubber closure of glass equipments used in aseptic techniques and numerous objects that can withstand high temperature and pressure

## Specification

Model	LTTA-A10	LTTA-A11	LTTA-A12	LTTA-A13
Capacity	20 L	24 L	35 L	50 L
Number of sterilizing trays	3	3	4	4
Sterilization Temperature	134 °C			
Sterilization Pressure	0.22 Mpa			
Timer Range	0-60 Minutes			
Heat Average	$\leq \pm 1^{\circ}\text{C}$			
Chamber Dimension (Φ×D)	250×420 mm	280×520 mm	300×500 mm	340×550 mm
Net Weight	38 Kgs	44 Kgs	80 Kgs	93 Kgs
Gross Weight	43 Kgs	50 Kgs	95 Kgs	112 Kgs
Power Supply	1.5 KW/AC220 V 50 Hz	1.5 KW/AC220 V 50 Hz	2.7 KW/AC220 V 50 Hz	2.7 KW/AC220 V 50 Hz