



Universal Pipette Tips UPT118L

Description :

Specifications :

Capacity	5000 µL
Length of tube	123.60 mm
Packaging	100 Tips/Bag, 10 Bags/Case
Quantity per Case	1000

Features :

Universal compatibility

Material: Polypropylene

Autoclavable at +121°C

DNase, RNase and Human DNA free

Micro Point Design

Pre Sterilized

Long reach design

Optional Color “ Yellow (UPT106L, UPT107L, UPT108L, UPT109L, UPT110L and UPT111L)

Optional Color “ Blue (UPT112L, UPT113L, UPT114L, UPT115L, UPT116L, UPT117L and UPT118L)



Labodam Equipment Ltd

18a Melton Road Leicester LE4 5EA
United Kingdom

www.labodam.com // info@labodam.com