



**VARIABLE VOLUME  
FULLY AUTOCLAVABLE  
PIPETTES LD-VVP110L**

## Description:

- Autoclavable at 121 °C.
- Material: Polycarbonate.
- Available in 0.1 µL to 5000 µL volume.
- Ergonomic design provides excellent operating experience.
- Easy-to-read volume display.
- Trouble-free calibration and maintenance.

## Features:

- Autoclavable at 121 °C
- Material: Polycarbonate
- Available in 0.1 µL to 10000 µL volume
- Ergonomic design provides excellent operating experience
- Easy-to-read volume display
- Trouble-free calibration and maintenance

## Specification:

Model no.	LD-VVP110L
Volume Range	2000 - 10000 µL
Increment	0.1 mL
Test Volume	2000 µL
	5000 µL
	10000 µL
Accuracy error	± 0.60 %
	± 1.20 %
	± 1.50 %
Precision error	± 0.20 %
	± 0.30 %
	± 0.40 %