



*Variable Volume Fully Autoclavable  
Pipettes VVP100L*

[www.labodam.com](http://www.labodam.com) // [info@labodam.com](mailto:info@labodam.com)

## Description :

## Specifications :

Volume Range	0.1 - 2.5 $\mu\text{L}$
Increment	0.05 $\mu\text{L}$
Test Volume	2.5 $\mu\text{L}$
	1.25 $\mu\text{L}$
	0.25 $\mu\text{L}$
Accuracy error	2.50 %
	3.00 %
	12.00 %
Precision error	2.00 %
	3.00 %
	6.00 %

## Features :

Autoclavable at 121  $^{\circ}\text{C}$

Material: Polycarbonate

Available in 0.1  $\mu\text{L}$  to 5000  $\mu\text{L}$  volume

Ergonomic design provides excellent operating experience

Easy-to-read volume display

Trouble-free calibration and maintenance



**Labodam Equipment Ltd**

**18a Melton Road Leicester LE4 5EA  
United Kingdom**

**[www.labodam.com](http://www.labodam.com) // [info@labodam.com](mailto:info@labodam.com)**