



*Variable Volume Fully Autoclavable  
Pipettes VVP102L*

[www.labodam.com](http://www.labodam.com) // [info@labodam.com](mailto:info@labodam.com)

## Description :

## Specifications :

Volume Range	2 - 20 µL
Increment	0.5 µL
Test Volume	20 µL
	10 µL
	2 µL
Accuracy error	0.90 %
	1.20 %
	3.00 %
Precision error	0.40 %
	1.00 %
	2.00 %

## Features :

Autoclavable at 121 °C

Material: Polycarbonate

Available in 0.1 µL to 5000 µL volume

Ergonomic design provides excellent operating experience

Easy-to-read volume display

Trouble-free calibration and maintenance



**Labodam Equipment Ltd**

18a Melton Road Leicester LE4 5EA  
United Kingdom

[www.labodam.com](http://www.labodam.com) // [info@labodam.com](mailto:info@labodam.com)