



*Variable Volume Fully Autoclavable  
Pipettes VVP108L*

## Description :

## Specifications :

Volume Range	200 - 1000 $\mu$ L
Increment	5 $\mu$ L
Test Volume	1000 $\mu$ L
	500 $\mu$ L
	200 $\mu$ L
Accuracy error	0.60 %
	0.70 %
	0.90 %
Precision error	0.20 %
	0.25 %
	0.30 %

## Features :

Autoclavable at 121  $\hat{A}$  $^{\circ}$ C

Material: Polycarbonate

Available in 0.1  $\hat{A}$  $\mu$ L to 5000  $\hat{A}$  $\mu$ L volume

Ergonomic design provides excellent operating experience

Easy-to-read volume display

Trouble-free calibration and maintenance



**Labodam Equipment Ltd**

18a Melton Road Leicester LE4 5EA

United Kingdom

[www.labodam.com](http://www.labodam.com) // [info@labodam.com](mailto:info@labodam.com)